



apex locators - micromotors



## Endy NT2

micromotor with light




Learn about the other products in our [ionyx](#) range

### Locapex Five

apex  
locator

Apex locator  
with instant  
measurement  
system



### Endy 6200

micromotor  
+ apex locator

Programmable  
endodontic  
micromotor with  
built-in apex locator



# Endy NT2

## micromotor with light

### rapid motor shutdown

As soon as the programmed torque is achieved, the motor is quickly shut down by means of an electronic lock, limiting the risk of breakage of the instrument.

### interruptible slow reverse

This 1-second slow-reverse technology used by IONYX since 1996, lets you remove the instrument without sending dental debris back towards the apex.

### long Li-ion battery life

The Endy NT2 lets you work on battery power for 1 hour with an average torque of 1.5 N.cm (40 root canal treatments on average reported).

### quick charger

The device's quick charger fully recharges your battery in 2 hours. Your battery is protected against short circuits or any heating during the charging process, ensuring complete safety. A red indicator light indicates that the battery is in safety lock-out.

### 10-memory setting

The Endy NT2 has a 10-memory setting for programming as many torques and speeds corresponding to each of the files used. Programming is easy and straightforward.

front view




side view




The Endy NT2 electronic casing has a 2-years warranty.  
The contra-angle and the battery have a 1-year warranty.

## Product strengths

a balanced compact motor  
with LED that can be autoclaved at 134°



The white LED under the contra-angle illuminates the operating field without casting a shadow or altering visibility.



a wide range of speeds  
with a single contra-angle

The Endy NT2 lets you set your file speed between 250 and 1200 revolutions per minute using the same contra-angle. A single machine can suffice for widening all root canals.

consistent torque and speed

The speed of the Endy NT2 remains consistent regardless of the force with which the file is used in the root canal. It never varies, even if the battery is less than 50% charged. The torque is always controlled in order to preserve your files.

*Your partner  
for over 20 years*

ionyx

34 rue de Fleurenne  
F-33290 BLANQUEFORT

T. +33 5 56 35 34 27  
F. +33 5 56 95 16 20

ionyx@ionyx.eu  
www.ionyx.eu

